




Introducing AppliedVR

The VR headset is a tool to enhance a patient's experience and make them as comfortable and relaxed as possible while at your center. Ask your patient if they would like to try the headset to experience a 3D video to help take their mind off their anxiety or pain and put them at ease.

Virtual Reality is a completely safe, non-invasive experience designed for all ages. Patients who have used it have said it was effective in helping them feel more relaxed. Patients can remove the headset at any time.

HOW TO ADMINISTER APPLIEDVR TREATMENT

- 1 Provide your patient with the information sheet that lists risks and benefits.
- 2 Turn on the mobile phone by holding down the button on the right side until it vibrates and powers on. Plug in the headphones at the base of the phone.
- 3 Click on the AppliedVR application on the home screen. Select an experience from the menu. Follow the prompt: "To open this application, insert your device into your Gear VR."
- 4 With the screen facing the headset, insert the phone in the headset using the **Accessory Interface Connector**. Press the phone into the headset until the lock clicks. Make sure the phone is securely fastened.



- 5 Before placing the headset on your patient's head, point out the **focus dial** located on the top of the headset and **volume button** on the right side to adjust as needed.
- 6 Ensure your patient has the hygienic hair cap on and place the headset on your patient's head. You can adjust the fit using the Velcro straps on either side and top of the headset. Then place the headphones on your patient.
- 7 Press the **back button**  on the right side of the headset to reorient the experience. The app may take a few seconds to load. If your patient changes positions, press the **back button** again to re-center their view.




Note: If the application does not load, remove the phone and reinsert the phone in the headset to ensure the phone is plugged in correctly.

- 8 After use, remove the phone from the headset by releasing the latch opposite the **Connector** to avoid draining the phone battery.

If patient reports motion sickness, claustrophobia or vertigo, please discontinue use of headset.

BETWEEN PATIENT USES

- 1 The Gear VR headset is fitted with the disposable hygienic liner. Remove the liner and dispose.
- 2 Place the darker, black side of a new liner on the Velcro part of the headset.



- 3 Completely wipe down the headset (plastic housing, lenses, straps, phone, headphones) using a Sani-wipe. Let sit for at least 2 minutes before next use.

Please contact AppliedVR at 844-857-0010 if you have any other questions or need assistance.